



DIMENSIONS IN MILLIMETRES DO NOT SCALE

Drawing & contents remain the property of Invenco Pty Ltd unless otherwise assigned under terms of engagement or other relevant agreement.

PARTS LIST				
ITEM	PART NUMBER	DESCRIPTION	QTY	MATERIAL
1		60ml Syringe - Barrel	1	HDPE
2		60ml Syringe - Cap	1	HDPE
3		60ml Syringe - Plunger	1	HDPE
4		60 Syringe - Ring	1	HDPE
5		60ml Syring - Piston	1	LDPE

Designed By Dave Scott	Checked By	Approved By	Date 10/02/2011	Client	Project
CPS Pty Ltd 52-54 Alexander Ave, Taren Point, NSW 2229, Australia Ph: 02 9531 1122 Fax: 02 9531 1144			Description 60ml Syringe		
			Part No.	Sheet 1 / 3	Rev.